



Impact of building material emissions on indoor air quality: simultaneous quantification of VOCs and formaldehyde in indoor air and at the air/material interface

Delphine Bourdin^{1,2}

Pierre Mocho³

Christophe Cantau¹

Valérie Desauziers²

¹ Nobatek – France

² L.G.E.I Mining School of Alès – France

³ LaTEP University of Pau - France

Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

VOCs and Formaldehyde in indoor air

French Legislation

Indoor air quality

- Decree n°2012-14: measurement of some indoor air pollutants in buildings open to the public (formaldehyde, benzene, CO₂)
- Decree n°2011-1727: guide-values for formaldehyde (10 µg/m³) and benzene (2 µg/m³)

Labeling of all building materials according to their VOCs emissions

- Decree n°2011-321: labeling of building products, flooring, wall covering, paints and varnishes on their volatile pollutants emissions (Formaldehyde, acetaldehyde, toluene, tetrachloroethylene, xylene, trimethylbenzene, dichlorobenzene, ethylbenzene, 2-butoxyethanol, styrene, TVOC)

Main source of indoor pollutants: building and decoration materials, furniture ¹



¹ Brown, Indoor Air 12 (2002) 55-63



Efficient tools to study both **indoor air** and **air/material interface** are needed

Indoor air
diagnosis

Identification of materials
responsible for
VOC pollution

Modeling of
Indoor air quality

Lyon, 26 septembre 2012



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Analytical strategies

Indoor air

On-site material emissions



DNPH
Cartridges




Adsorbent tubes




Radiello®



FLEC®
(ISO 16000-10)



Shinohara et al., 2009
Blondel and Plaisance, 2010



VOCs (GC)



Aldehydes (HPLC)

Solid-Phase Micro Extraction (SPME)



Desorption directly in a GC injector

Development of an analytical method for
the simultaneous quantification of VOCs
and aldehydes

Lyon, 26 septembre 2012



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Analytical Performances

Indoor air

On-site material emissions



Instantaneous sampling

SPME fiber

Desauziers, Auguin,
Techniques de
l'ingénieur, 2012

✓ Mix of 9 compounds:

- Formaldehyde*
- Acetaldehyde*
- Toluene*
- p-xylene*
- Styrene*
- Alpha-pinene°
- Dichlorobenzene*
- Tetrachloroethylene*
- Hexanal°

- * Concerned by material labeling
- ° Wooden construction

Patent
pending...

SPME – emission cell coupling

Non-equilibrium → First Fick's law

$$n_f = U \cdot (C_i \cdot t)$$

n_f : extracted amount on the fiber

U : uptake rate

C_i : concentration in air ($\mu\text{g} \cdot \text{m}^{-3}$)

t : extraction time (min)




Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Lyon, 26 septembre 2012

Analytical Performances: calibration curves (FID detection)



Formaldehyde

- ✓ Good linearity for the 9 compounds
- ✓ Should be possible to analyse higher concentrations for other applications...

Compound	R ²	Linearity range for a 5 minute SPME extraction ($\mu\text{g.m}^{-3}$)
Acetaldehyde	0.9178	528
Formaldehyde	0.9959	480
Toluene	0.9720	920
Tetrachloroethylene	0.9688	1660
P-xylene	0.9858	1480
Styrene	0.9746	1460
Alpha-pinene	0.9653	1900
hexanal	0.9699	1500
1,2-dichlorobenzene	0.9616	2040
AVERAGE		1328



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Lyon, 26 septembre 2012

Analytical Performances: Limits of detection (LOD), limits of quantification (LOQ) and reproducibility

SPME: 20 minute exposure

Concentration of each compound in the standard gas: $2.4 \mu\text{mol.m}^{-3}$

Results in $\mu\text{g/m}^3$

Compounds	GC/FID			GC/MS – SIM mode	
	% RSD (n = 6)	LOD	LOQ	LOD (S/N = 3)	LOQ (S/N = 10)
Toluene	17%	3.6	11.8	0.012	0.042
Tetrachloroethylene	14%	4.7	15.8	0.014	0.048
P-xylene	8%	2.1	6.8	0.027	0.091
Styrene	5%	1.4	4.7	0.013	0.043
Alpha-pinene	11%	3.3	10.9	0.049	0.164
1,2-dichlorobenzene	9%	2.8	9.5	0.011	0.037
Formaldehyde	16%	4.2	14.0	0.005	0.016
Acetaldehyde	22%	11.3	37.8	0.051	0.170
Hexanal	13%	3.7	12.3	0.124	0.414
AVERAGE		4.1	13.7	0.034	0.113




Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Lyon, 26 septembre 2012

Influence of air relative humidity




Formaldehyde



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71


Alpha-pinene



Lyon, 26 septembre 2012

Comparison with the standard method

	NF ISO 16000-3 (2002) DNPH Cartridges	SPME
Sampling	3.5 hours for indoor atmospheres 15 minutes for material emissions	Instantaneous sampling 15 minute SPME extraction
Analyses	HPLC/UV	GC/FID



Material emissions in environmental chamber



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Indoor air



Lyon, 26 septembre 2012

Case Study



Siège Social

NOBATEK

67, rue de Mirambeau

64600 ANGLET

Tél. : 05 59 03 61 29

Fax. : 05 59 63 55 41

www.nobatek.com

Site Ecocampus

ENSAM

Esplanade des Arts et Métiers

33405 TALENCE

Tél. : 05 56 84 63 70

Fax. : 05 56 84 63 71

Office building located in Mont-de-Marsan (40)
Occupied since one year
Built in a HQE approach

Lyon, 26 septembre 2012

Case Study



- ✓ Meeting room of 40 m²
- ✓ Study of indoor air by vacuum vial
- ✓ Study of surface's concentration by emission cell
 - Floor
 - Tables
 - Chairs
 - Doors
 - Ceiling
 - Roller blind



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

➡ Identification of terpenes, aldehydes and alkanes

➡ 4 compounds belonging to our model compounds were quantified:
alpha-pinene, formaldehyde, acetaldehyde and hexanal

Lyon, 26 septembre 2012

Formaldehyde


Material emissions

Surface	Air interface concentration ($\mu\text{g.m}^{-3}$)
Table	42.8
Floor	15.3
Chair	6.7
Ceiling (painted part)	7.7
Wall (paint)	10.1
Door	16.5
Ceiling (ventilation area)	11.8
Roller blind	10.2

Indoor concentration
($\mu\text{g.m}^{-3}$)

7.0

Measurement of air exchange rates by CO_2 injection



Indoor air quality
modeling



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Lyon, 26 septembre 2012

Indoor Air Modeling

Mass Balance

$$\frac{dC_i}{dt} = \sum_{j=1}^m Q_{ji} + \lambda_{air} f C_{ext} - \lambda_{air} C_i - U_{dep} \left(\frac{A}{V} \right) C_i$$

Indoor Concentration

Material Input

Outdoor Concentration

Room Ventilation

Sink Effect

C_i : indoor air concentration (ppb)

t : time (s)

Q_{ji} : emission rate (ppb.s⁻¹)

λ_{air} : air exchange rate (s⁻¹)

C_{ext} : outdoor concentration (ppb)

U_{dep} : deposit velocity (m.s⁻¹)

A : surface (m²)

V : room volume (m³)

f : penetration rate (0 < f < 1)

$$C_i(t) = C_0 e^{-(U_{dep} \frac{A}{V} + \lambda_{air})t} + \frac{\sum_{j=1}^m Q_{ji} + \lambda_{air} C_{ext}}{U_{dep} \frac{A}{V} + \lambda_{air}} \left(1 - e^{-(U_{dep} \frac{A}{V} + \lambda_{air})t} \right)$$



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71


Lyon, 26 septembre 2012

Indoor air modeling: application to formaldehyde

$$U_{\text{dep}} = 0 \text{ m.s}^{-1}$$
$$C_{\text{ext}} = 1 \text{ }\mu\text{g.m}^{-3}$$
$$\lambda_{\text{air}} = 0.98 \text{ vol/h}$$

Calculated formaldehyde concentration: $6.4 \text{ }\mu\text{g.m}^{-3}$

Measured formaldehyde concentration: $7.0 \pm 1.0 \text{ }\mu\text{g.m}^{-3}$



=> In accordance with the future French regulation ($10 \text{ }\mu\text{g.m}^{-3}$)



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

- Emission rate calculated from the air/material interface concentration
- Modeling considering the air exchange rate and the “real” material contribution
- Possibility to try different scenarii in order to evaluate the impact of various parameters of building management: ventilation, materials....

Lyon, 26 septembre 2012

Conclusion

- ✓ Fast, sensitive, easy to implement analytical method for the evaluation of **both VOCs and formaldehyde** was developed. Further comparisons with the standard methods will now be done.
- ✓ On-site determination of material emissions
- ✓ Modeling of indoor air quality:
 - First experiment in a new building very encouraging
 - The box model will be adjusted with a 6 month study in a newly built high school
 - Predictive tool of indoor air quality



Siège Social
NOBATEK
67, rue de Mirambeau
64600 ANGLET
Tél. : 05 59 03 61 29
Fax. : 05 59 63 55 41
www.nobatek.com

Site Ecocampus
ENSAM
Esplanade des Arts et Métiers
33405 TALENCE
Tél. : 05 56 84 63 70
Fax. : 05 56 84 63 71

Thank you for your attention

Questions?



Lyon, 26 septembre 2012